

PCR analysis of lambda phage

Judelson 9-97

1. Using a sterile pasteur pipette, carefully remove the desired plaque(s) from the plate. GENTLY expel the plaque into a 1.5 ml tube containing 200 μ l of 1/2 strength SM. Add 20 μ l chloroform and gently hit the tube with your fingers. Let sit 2 hr at room temp.
2. Remove 5 μ l of the phage suspension and place in PCR tube containing 10 μ l water. Heat 5 minutes at 98^oC.
3. Add 15 μ l of 2X concentrated PCR reaction mix.

Recipe for 2X mix (1 reaction)

3 μ l 10X PCR buifer

0.3 μ l 100 x BSA stock (10 mg/ml)

2.4 μ l 1.25 mM dNTP solution

0.24 μ l Taq polymerase

50 ng of each primer (0.05 μ l of 1 mg/ml stock).

4. Do amplification with appropriate primer set. For the GT10 primers: use 50C anneaiing, 30-35 cycles.