

PLAQUE LIFTS

1. Plate the phage using standard procedures. Aim for a maximum of 40,000 pfu/plate for 150 mm plate, or 12,000 per 100 mm plate. The object is to obtain plaques that are barely touching each other. Use >1-day old plates (to avoid condensation smearing plaques and to obtain better adherence between top and bottom agar).
2. Incubate at **37C** for ca. 8 hrs.
3. Chill the plates for 2 hr at 4C to firm up the top agar.
4. Transfer the phage to filters: Gently lay a nitrocellulose membrane on the plate for 2 minutes. Orient with ink.

If making duplicate lifts, let the second filter transfer for 5 mm. Wrap the plates and store at 4⁰C.

5. Denature the membrane by placing (plaque side up) on 1.5 M NaCl, 0.5 M NaOH, 2 min.
6. Neutralize by placing on 1.5 M NaCl, 0.5 M Tris pH 8.0, 5 min.
7. Rinse by submerging in 0.2 M Tris-Cl PH 7.5, 2 X SSC (or 3X SSPE) for 30 sec maximum.
8. Blot briefly on Whatman paper.
9. Oven bake 80C, 1.5 hr.
10. Hybridize as usual.
11. Picking positives:

Using a sterile pasteur pipette, carefully remove the plaque(s) from the original plate. GENTLY expel the plaque into a 1.5 ml tube containing 300 μ l SM. Add 20 μ l chloroform and gently hit the tube with your fingers. Let sit 6-8 hr at 4C (or 2 hr at room temp). Assaying 10⁷ pfu/plaque, plate dilutions and do secondary (and tertiary) screens. It will usually take 3 rounds to purify phage due to the ease by which they diffuse. Just in case, store the primary plaque suspension at 4C (NEVER freeze lambda phage suspensions).