

Plating Lambda Libraries

1. Preparation of host bacteria
 - a. Inoculate LB (+ 10mM MgSO₄ and 0.2% maltose added after autoclaving) with a single colony. Grow at 37C, shaking, for 4-6 hours (OD₆₀₀ max=1.0) or 30C overnight.
 - b. Pellet bacteria at full speed (setting 4) in table top centrifuge, 10 min.
 - c. Gently resuspend in 10 mM MgSO₄ at about OD₆₀₀=0.5.
2. Plating phage.

(scale up these volumes as needed for your experiment)

 - a. assuming that a single plaque will contain about 10⁵ phage (or a plate lysate 10⁹/ml), calculate the required dilutions for your experiment.
 - b. To a sterile 16 ml plastic tube, add the lambda phage (in 100 µl SM) to 200 µl bacteria (in 10 mM MgSO₄). Be sure not to transfer any chloroform to the mixture! Incubate at 37C for 15-20 min to allow the phage to attach to the cells.
 - c. Add 3 ml of top agar (NZY broth gives bigger and clearer plaques), melted and cooled to 48C (±2C). QUICKLY mix by inverting 4X and pour IMMEDIATELY on prewarmed bottom agar plates (NZY).
 - d. Wait a few minutes for the plates to solidify. Incubate the plates 7-12 hrs at 37C, inverted.

SM (1 liter)

5.8 g NaCl

2 g MgSO₄·7H₂O

6.05 g Tris base

water

pH to 7.5 with NaOH

Add 0.1 g gelatin

autoclave

NZY broth (1 liter)

5 g NaCl

2 g MgSO₄·7H₂O

10 g NZ amine (casein hydrolysate) water

pH to 7.5 with NaOH

for top agar: add agarose to 0.7%.

for bottom agar: add agar to 1.5% (standard applications) OR agarose (for DNA preps).